


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/28/17
Time Processed : 14:01:28
Injection Volume : 0.2
Sample Name : Pancreatic Polypeptide, human
Cat. No. : RP10398
Lot. No. : PE2511712
Theoretical MW : 4181.70
Observed MW : 4182.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH