


# Mass Spectrum



## Sample Information

Month-Day Processed : 10/23/20  
Time Processed : 20:24:06  
Injection Volume : 0.3  
Sample Name : Angiotensin II (3-8), human  
Cat. No. : RP11253  
Lot. No. : E8717720J  
Theoretical MW : 774.90  
Observed MW : 774.3

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H2O/50% MeOH