


# Mass Spectrum



## Sample Information

Acquired by: Gary  
Month-Day Processed: 06/10/19  
Time Processed: 12:06:18 PM  
Injection Volume: 0.3  
Sample Name: Bim BH3, Peptide IV  
Cat. No.: RP20339  
Lot. No.: PE2141906  
Theoretical MW: 3269.70  
Observed MW: 3269.5

Interface: ESI  
Nebulizing Gas Flow: 1.5L/min  
CDL Temp: 250  
Block Temp: 200

Equipment: GK11010007  
Interface Bias: +4.5 kV  
Drying Gas Flow: 5 L/min  
T.Flow: 0.2 ml/min  
B.conc: 50% H2O/50% MeOH