

Sample Name :Angiotensin II,human  
 Cat No.: RP11257  
 Lot No.: PE2251712

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
31.00	Pumps	Pump B Conc.	95	
31.01	Pumps	Pump B Conc.	5	
40.00	Pumps	Pump B Conc.	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Alltima™ C18 4.6 x 250 mm

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	10.734	9908	1548	0.080
2	10.890	37517	7443	0.303
3	11.381	32613	5270	0.264
4	11.727	9523	1806	0.077
5	12.197	14663	2064	0.118
6	12.546	30102	4937	0.243
7	13.057	5595	1585	0.045
8	13.148	23007	4446	0.186
9	13.342	4586	1693	0.037
10	13.436	67293	6810	0.544
11	13.667	22288	4460	0.180
12	13.753	24971	5118	0.202
13	13.874	63588	13936	0.514
14	14.026	11793451	1984623	95.298
15	14.217	34808	11423	0.281
16	14.317	22518	4574	0.182
17	14.474	41472	5937	0.335
18	14.600	48573	9074	0.392
19	14.737	32386	5863	0.262
20	14.979	20128	4134	0.163
21	15.058	14203	3442	0.115
22	15.129	10812	2857	0.087
23	16.295	11353	2445	0.092

Peak#	Ret. Time	Area	Height	Area %
Total		12375358	2095491	100.000