


Mass Spectrum



Sample Information
Month-Day Processed : 05/28/21
Time Processed : 7:24:47
Injection Volume : 0.4
Sample Name : Pep-1 (uncapped)
Cat. No. : RP10130
Lot. No. : E4402621E
Theoretical MW : 2848.30
Observed MW : 2848.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH