


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	03/20/21	:ESI		: GK11010007	
Time Processed :	14:57:17	Nebulizing Gas Flow :1.5L/min		Interface Bias	
Injection Volume :	0.5	CDL Temp :250		Drying Gas Flow	
Sample Name :	Fibronectin Active Fragment Control	Block Temp :200		T.Flow	
Cat No. :	RP20347			B.conc	
Lot No. :	E8090721C				
Theoretical MW :	504.50				
Observed MW :	504.2				