


# Mass Spectrum



## Sample Information

Acquired by : Gary  
 Month-Day Processed : 11/09/17  
 Time Processed : 14:48:32  
 Injection Volume : 0.2  
 Sample Name :  $\beta$ -Amyloid (1-40)  
 Cat No.: RP10004  
 Lot No.: PE6771710  
 Theoretical MW : 4329.82  
 Observed MW : 4330.40

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250°C  
 Block Temp : 200°C

Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50%H2O/50%MeOH