


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 07/20/18  
Time Processed : 17:47:20  
Injection Volume : 0.2  
Sample Name : Peptide YY (PYY) (3-36), human  
Cat. No. : RP10354  
Lot. No. : PE2861807  
Theoretical MW : 4049.48  
Observed MW : 4049.70

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH