


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 02/27/19
Time Processed : 07:51:46 AM
Injection Volume : 0.2
Sample Name : BNP (1-32), rat
Cat No. : RP11121
Lot No. : PE5331902
Theoretical MW : 3452.94
Observed MW : 3452.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH