


Mass Spectrum



Sample Information

Month-Day Processed : 06/16/21
Time Processed : 14:34:49
Injection Volume : 0.2
Sample Name : CEF3
Cat. No. : RP11765
Lot. No. : E4020421F
Theoretical MW : 911.10
Observed MW : 911.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH