


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/19/19
Time Processed : 18:00:42
Injection Volume : 0.3
Sample Name : HPV16 E7(86-93)
Cat. No. : RP20499
Lot. No. : E3463319L
Theoretical MW : 815.10
Observed MW : 814.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH