


Mass Spectrum



Sample Information		Equipment	
Month-Day Processed :	08/13/20	Interface	: ESI
Time Processed :	7:18:17	Interface Bias	: +4.5kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	Endomorphin-1	T.Flow	: 0.2ml/min
Cat No. :	RP10925	B.conc	: 50%H2O/50%MeOH
Lot No. :	E6240220H		
Theoretical MW :	610.70		
Observed MW :	610.2		