


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	03/05/19	Interface Bias	: +4.5kV
Time Processed :	11:24:30 AM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.2	T.Flow	: 0.2ml/min
Sample Name :	TFLLR	B.conc	: 50% $\text{H}_2\text{O}$ /50% $\text{MeOH}$
Cat No. :	RP20354		
Lot No. :	PE8851902		
Theoretical MW :	647.81		
Observed MW :	647.8		