


Mass Spectrum



Sample Information
Month-Day Processed : 06/18/21
Time Processed : 18:01:29
Injection Volume : 0.4
Sample Name : Melanotan (MT)-II
Cat. No. : RP17323
Lot. No. : E9420021D
Theoretical MW : 1024.18
Observed MW : 1023.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH