


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 05/16/18  
Time Processed : 18:48:02  
Injection Volume : 0.2  
Sample Name : [Arg<sup>8</sup>]-Vasotocin  
Cat. No. : RP10061  
Lot. No. : PE5011804  
Theoretical MW : 1050.22  
Observed MW : 1049.90

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH