


Mass Spectrum



Sample Information
Month-Day Processed : 01/15/21
Time Processed : 8:29:34 AM
Injection Volume : 0.3
Sample Name : LL-37, scrambled
Cat No. : RP20332
Lot No. : E7107421A
Theoretical MW : 4493.30
Observed MW : 4492.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH