


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 08/31/18
Time Processed : 22:34:46
Injection Volume : 0.2
Sample Name : Dynorphin A(1-13),porcine
Cat. No. : RP20492
Lot. No. : PE3811808
Theoretical MW : 1604.00
Observed MW : 1604.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH