


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 07/20/17  
Time Processed : 15:39:04  
Injection Volume : 0.2  
Sample Name : Angiotensin II (3-8), human  
Cat No. : RP11253  
Lot No. : PE9901707  
Theoretical MW : 774.90  
Observed MW : 775.00

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250 °C  
Block Temp : 200 °C

## Interface Bias

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH