


Mass Spectrum



Sample Information
Month-Day Processed : 01/30/21
Time Processed : 2:42:41 PM
Injection Volume : 0.5
Sample Name : Phe-Met-Arg-Phe, amide
Cat No. : RP10818
Lot No. : E4424421A
Theoretical MW : 598.76
Observed MW : 598.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H₂O/50%MeOH