


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 01/08/17
Time Processed : 13:45:26
Injection Volume : 1
Sample Name : Insulin alpha-chain (1-13)
Sample ID : U7815CA050-1
Theoretical MW : 1640.00
Observed MW : 1640.20

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH