


Mass Spectrum



Sample Information
User : Gary
Month-Day Processed : 03/26/18
Time Processed : 10:16:16
Injection Volume : 0.2
Sample Name : Gastrin-1, human
Cat No. : RP12740
Lot No. : PE9561803
Theoretical MW : 2098.22
Observed MW : 2098.10

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250 °C
Block Temp : 200 °C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH