


Mass Spectrum



Sample Information

Month-Day Processed : 02/07/21
Time Processed : 1:40:50 AM
Injection Volume : 0.3
Sample Name : Arg-Gly-Asp-Ser
Cat. No. : RP10861
Lot. No. : E3766321B
Theoretical MW : 433.42
Observed MW : 433.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH