


Mass Spectrum



Sample Information
Month-Day Processed : 07/13/21
Time Processed : 11:10:37
Injection Volume : 0.2
Sample Name : C-Peptide, dog
Cat No. : RP11115
Lot No. : E1542021G
Theoretical MW : 3174.49
Observed MW : 3174.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH