


# Mass Spectrum



## Sample Information

Month-Day Processed : 09/08/20  
 Time Processed : 02:49:10 PM  
 Injection Volume : 0.8  
 Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat  
 Cat No. : RP11121  
 Lot No. : E2516420H  
 Theoretical MW : 3452.94  
 Observed MW : 3452.0

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

Equipment : GK11010007  
 Interface Bias : +4.5kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2ml/min  
 B.conc : 50%H2O/50%MeOH