


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/05/19
Time Processed : 19:39:56
Injection Volume : 0.3
Sample Name : Peptide YY (PYY) (3-36), human
Cat No. : RP10354
Lot No. : E9593119I
Theoretical MW : 4049.48
Observed MW : 4048.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH