


Mass Spectrum



Sample Information
Month-Day Processed : 09/26/20
Time Processed : 07:38:23 AM
Injection Volume : 0.3
Sample Name : PACAP (1-27), human, ovine, rat
Cat No. : RP10334
Lot No. : E7997820I
Theoretical MW : 3147.62
Observed MW : 3146.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH