


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 08/02/20
Time Processed : 7:28:57
Injection Volume : 0.4
Sample Name : Fibronectin CSI Peptide
Cat No. : RP10838
Lot No. : E1435420G
Theoretical MW : 872.96
Observed MW : 872.3

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : -4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH