


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/22/16
Time Processed : 21:37:54
Injection Volume : 1
Sample Name : Melanotan (MT)-II
Cat No. : RP17323
Lot No. : PE3351612
Theoretical MW : 1024.18
Observed MW : 1023.55

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH