


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/15/20
Time Processed : 12:45:27
Injection Volume : 0.3
Sample Name : Bim BH3, Peptide IV
Cat No. : RP20339
Lot No. : E2881220B
Theoretical MW : 3269.70
Observed MW : 3269.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH