


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/14/17
Time Processed : 21:58:10
Injection Volume : 0.2
Sample Name : Calcineurin Substrate
Cat No. : RPI0139
Lot No. : PE3331710
Theoretical MW : 2112.40
Observed MW : 2112.15

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH