


Mass Spectrum



Sample Information

User : Gary
 Month-Day Processed : 04/02/19
 Time Processed : 10:17:30 PM
 Injection Volume : 0.2
 Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat
 Cat No. : RP11121
 Lot No. : PE7791904
 Theoretical MW : 3452.94
 Observed MW : 3453.0

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment : GK11010007
 Interface Bias : +4.5kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2ml/min
 B.conc : 50%H2O/50%MeOH