


Mass Spectrum



Sample Information
Month-Day Processed : 06/04/21
Time Processed : 10:47:40
Injection Volume : 0.2
Sample Name : Dynorphin A(1-13),porcine
Cat. No. : RP20492
Lot. No. : E3961721E
Theoretical MW : 1604.00
Observed MW : 1604.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH