


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/13/19
Time Processed : 06:54:02 PM
Injection Volume : 0.3
Sample Name : PACAP (1-38), human, ovine, rat
Cat No. : RP10335
Lot No. : E4707819I
Theoretical MW : 4534.26
Observed MW : 4533.6

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH