


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	01/18/20	Interface Bias	: +4.5kV
Time Processed :	20:56:23	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Kemptide	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP10287		
Lot No. :	E1559220A		
Theoretical MW :	771.90		
Observed MW :	771.6		

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200