


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	06/19/19	Interface Bias	: +4.5kV
Time Processed :	09:55:50 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.4	T.Flow	: 0.2ml/min
Sample Name :	[Glu1]-Fibrinopeptide B	B.conc	: 50%H2O/50%MeOH
Cat No. :	RP20439		
Lot No. :	E3923119F		
Theoretical MW :	1570.60		
Observed MW :	1570.4		