


Mass Spectrum



Sample Information
Month-Day Processed : 01/05/21
Time Processed : 10:04:27 AM
Injection Volume : 0.3
Sample Name : InsB 9-23
Cat No. : RP20391
Lot No. : E7173820L
Theoretical MW : 1645.90
Observed MW : 1645.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH