


# Mass Spectrum



## Sample Information

Month-Day Processed : 09/28/20  
Time Processed : 20:44:17  
Injection Volume : 0.3  
Sample Name : Mastoparan  
Cat No. : RP10465  
Lot No. : E4901520I  
Theoretical MW : 1478.80  
Observed MW : 1478.7

## Interface

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

## Equipment

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH