


Mass Spectrum



Sample Information

Month-Day Processed :	04/30/21	Equipment	: GK11010007
Time Processed :	7:32:13 AM	Interface Bias	: +4.5 kV
Injection Volume :	0.7	Drying Gas Flow	: 5 L/min
Sample Name :	[Glu1]-Fibrinopeptide B	T.Flow	: 0.2 ml/min
Cat. No. :	RP20439	B.conc	: 50% H2O/50% MeOH
Lot. No. :	E7689821D		
Theoretical MW :	1570.60		
Observed MW :	1570.4		