


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 12/06/19  
Time Processed : 21:21:10  
Injection Volume : 0.4  
Sample Name : Melanotan (MT)-II  
Cat No. : RP17323  
Lot No. : E0098419L  
Theoretical MW : 1024.18  
Observed MW : 1024.2

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH