


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 01/21/19
Time Processed : 08:27:03 PM
Injection Volume : 0.2
Sample Name : 43Gap 26.Connexin Mimetic
Cat No. : RP20274
Lot No. : PE6141901
Theoretical MW : 1550.79
Observed MW : 1551.0

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH