


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 11/22/19  
Time Processed : 08:28:16 PM  
Injection Volume : 0.4  
Sample Name : Somatostatin  
Cat. No. : RP10226  
Lot. No. : E3616419K  
Theoretical MW : 1637.88  
Observed MW : 1637.7

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50% H<sub>2</sub>O/50% MeOH