


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 09/18/19
Time Processed : 03:18:03 PM
Injection Volume : 0.3
Sample Name : HPV16-E711-20 epitope
Cat No. : RP20285
Lot No. : E4869619I
Theoretical MW : 1210.38
Observed MW : 1209.6

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH