


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 12/09/17
Time Processed : 11:24:43
Injection Volume : 5
Sample Name : BNP (1-32), rat
Cat. No. : RP11121
Lot. No. : PE2031711
Theoretical MW : 3452.94
Observed MW : 3453.00

Interface

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias

Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH