


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 03/06/19
Time Processed : 10:17:22 AM
Injection Volume : 0.2
Sample Name : Brain Natriuretic Peptide (BNP) (1-32), rat
Cat No. : RP11121
Lot No. : PE7951902
Theoretical MW : 3452.94
Observed MW : 3452.5

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H₂O/50%MeOH