


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 02/01/20
Time Processed : 07:52:04 PM
Injection Volume : 0.3
Sample Name : Melanin Concentrating Hormone, salmon
Cat No. : RP10662
Lot No. : E7232620B
Theoretical MW : 2099.48
Observed MW : 2098.8

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH