


Mass Spectrum



Sample Information

|                       |                 |                 |                  |
|-----------------------|-----------------|-----------------|------------------|
| User :                | Gary            | Equipment       | : GK11010007     |
| Month-Day Processed : | 07/24/19        | Interface Bias  | : +4.5kV         |
| Time Processed :      | 02:58:34 PM     | Drying Gas Flow | : 5 L/min        |
| Injection Volume :    | 0.3             | T.Flow          | : 0.2ml/min      |
| Sample Name :         | Arg-Gly-Asp-Cys | B.conc          | : 50%H2O/50%MeOH |
| Cat No. :             | RP10860         |                 |                  |
| Lot No. :             | E9047019G       |                 |                  |
| Theoretical MW :      | 449.48          |                 |                  |
| Observed MW :         | 449.2           |                 |                  |