


# Mass Spectrum



Sample Information

Month-Day Processed : 10/21/20  
 Time Processed : 04:53:46 PM  
 Injection Volume : 0.3

Sample Name : Peptide YY (PYY) (3-36), human  
 Cat No. : RP10354  
 Lot No. : E1765220J

Theoretical MW : 4049.48  
 Observed MW : 4048.8

Interface : ESI  
 Nebulizing Gas Flow : 1.5L/min  
 CDL Temp : 250  
 Block Temp : 200

Equipment : GK11010007  
 Interface Bias : +4.5 kV  
 Drying Gas Flow : 5 L/min  
 T.Flow : 0.2 ml/min  
 B.conc : 50% H2O/50% MeOH