


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 01/15/17  
Time Processed : 13:01:47  
Injection Volume : 1  
Sample Name : Insulin alpha-chain (1-13)  
Cat No. : RP20501  
Lot No. : PE8301701  
Theoretical MW : 1640.00  
Observed MW : 1639.90

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250°C  
Block Temp : 200°C

Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H2O/50%MeOH