


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 11/30/17
Time Processed : 16:39:07
Injection Volume : 0.2
Sample Name : [Met⁵]-Enkephalin
Cat. No.: RP10886
Lot. No.: PE2091711
Theoretical MW : 573.66
Observed MW : 573.20

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH