


Version 1.1

Update: 9/28/2019

# Certificate of Analysis

## Product Information

<b>Product Name:</b>	Human CD4 Antibody (M5A8) [PE], mAb, Mouse
<b>Cat. No.:</b>	V02702-25T
<b>Quantity:</b>	25 tests
<b>Presentation:</b>	PBS, pH 7.4, containing gelatin and 0.03% ProClin300.
<b>Storage:</b>	The antibody should be stored at 4°C and protected from prolonged exposure to light. <b>Do not freeze.</b>
<b>Purification:</b>	The antibody was purified by affinity chromatography and conjugated with PE under optimal conditions. The solution is free of unconjugated PE.
<b>Lot No.:</b>	PP200601
<b>Concentration:</b>	40 µg/mL
<b>Date of Manufacturing:</b>	06/29/2020
<b>Expiry Date:</b>	06/29/2021
<b>Quality control:</b>	Flow cytometric analysis of antibody surface-stained cells.



Human peripheral blood lymphocytes were stained with or without Human CD4 Antibody (M5A8) [PE].

**Date:** 06/30/2020

**Approved by:** *Xianli Yang*

**Quality Department**

*For laboratory research use only. Direct human use, including taking orally and injection and clinical use are forbidden.*